

# Telehealth Mobile App

APPLE INSTALLATION GUIDE | MAY 2020



**CONNECT***withUs.*

**eVerO Corporation – Corporate Headquarters**

48 South Service Road

Suite 200

Melville, NY 11747 · USA

1.866.GO.EVERO (1.866.463.8376)

 <http://www.evero.com>


 [support@evero.net](mailto:support@evero.net)

Copyright © 2020 eVerO Corporation. All rights reserved.

digitalAGENCY™ logo is trademarked, and eVerO is a registered trademark of eVerO Corporation. All other trademarks used herein are the properties of their respective owners.


This installation guide will provide you with step-by-step instructions on how to install Telehealth Mobile to your personal device. Since our app is not available through the App Store (iOS) or Google Play, additional steps are required for installation.

Using your Apple device, you will need to install the app from our in-house app install system. To do this, Navigate to the eVero website (*located at eVero.com*) and click on **'Our Apps'** located in the navigation bar to the top of the screen.



Once you have accessed the **'Our Apps'** page, you will find a brief overview of the apps eVero has to offer. To access further information in regard to downloading the Telehealth Mobile App click on the designated button.

On this page, you will now have access to the direct download link for the app by clicking on the appropriate download button.




When prompted, click **Install** to begin.

To complete the additional steps required for Apple installation, tap the Settings icon that is typically found on the Home screen.




Once you have opened the **Settings** screen, select **General** and tap **Device Management**.

You will see **Evero Corporation** listed under **Enterprise App**. Tap on **Evero Corporation**.



Click on **Trust "Evero Corporation"** as shown below. This will allow your mobile device to install Telehealth Mobile through a source outside of the App Store.



Copyright @ 2020 eVero Corporation. All rights reserved.

digitalAGENCY™ logo is trademarked, and eVero is a registered trademark of eVero Corporation. All other trademarks used herein are the properties of their respective owners.